

Preparation of Angiotensin and Vasoactive Peptide Stocks

Materials:

- Angiotensin I human acetate salt (Sigma A9650-1MG)
- Vasoactive intestinal peptide fragment 1-12 human, porcine, rat (Anaspec AS-24217, 1 mg)
- Mobile phase A: 0.1% formic acid in water (Optima LC/MS grade)

1. Make 1 mM stock solutions of both peptides in original container:
 - 1 mM Angiotensin Stock: Dissolve 1 mg solid in 772 μL Solvent A
 - 1 mM Vasoactive Stock: Dissolve 1 mg solid in 702 μL Solvent A
2. From the 1 mM stock solutions, make a 1 pmol/ μL stock solution containing both peptides.
 - In Lo-Bind tube, to **998 μL** of Solvent A:
 - (1) Add 1 μL of 1 mM Angiotensin stock and
 - (2) Add 1 μL of 1 mM Vasoactive stock.
 - Total = 1000 μL of 1 pmol/ μL AV Standard stock solution
3. Aliquot stock solution into smaller volumes (25 μL). Each lab member will keep their own stock of concentrated 1 pmol/ μL AV standards and use this to dilute AV standards to concentrations for spiking into proteomic samples/mass spec standards.
4. For mass spec standards, make 10x dilution from the concentrated 1 pmol/ μL stock to make 100 fmol/ μL .

Revised: 07/30/2022