

Prep Control (to ensure sample quality):

HEK293T Sol Proteome

***In Vitro* HHS-465 100uM 1hr rt treatment**

SILAC: Probe/Probe

For detection of active-site peptides.

***when measuring inhibitor dry stocks for DMSO solutions, do not try to measure less than 2 mg of compound. The error from measuring small amounts of compound can drastically affect the final concentrations tested in your assay.**

A-1. START HERE if you are performing inhibitor competition

1. Thaw 1 mg proteome aliquots on ice for sample and control
 - *proteomes are in PBS with no protease or phosphatase inhibitors
 - *proteome concentrations are 2.3 mg/mL to achieve 1 mg using volumes specified below
2. Use 427 μ L of proteome from #1 for the preparation below.
3. Add 5 μ L of 100x compound (1 mM solution in DMSO), and incubate 30 min, 37 °C
4. Follow 3 and 4 in A-2

A-2. START HERE if you are only performing probe labeling

1. Thaw 1mg proteome aliquots on ice for sample and control
 - *proteomes are in PBS with no protease or phosphatase inhibitors
 - *proteome concentrations are 2.3 mg/mL to achieve 1 mg using volumes specified below
2. Use 432 μ L of proteome from #1 for the preparation below.
3. Add 5 μ L of 100x probe (or DMSO to control), gently mix tube.
4. Incubate 1 hr, room temperature (RT)
 - *The time may be adjusted based on the probe.

B. Click chemistry (use only with alkyne probe; for biotin probes go to procedure C)

***if you are making a master click chemistry mix (i.e. all compounds premixed before adding to individual samples), do not add CuSO₄ until just before you are ready to add master mix to samples.**

1. Create master-mix by adding:
 - a. Tag: Add 10 μ L of 50x (10 mM, 2 mg/410 μ L in DMSO) Biotin-N₃ or desthiobiotin-N₃, giving **200 μ M final**.
 - b. TCEP: Add 10 μ L of 50x (50 mM, 7.17 mg/500 μ L in H₂O) TCEP (**make fresh**), giving **1 mM final**, vortex.
 - c. Ligand: Add 33 μ L of 16.7x (1.67 mM) TBTA Ligand, giving **100 μ M final**, vortex **1.67 mM = 20 μ L 50x + 180 μ L DMSO + 800 μ L t-BuOH 8.55mg into 200 μ L in DMSO gives 50x stock**
 - d. CuSO₄: Add 10 μ L of 50x (50 mM, 12 mg in 960 μ L H₂O) CuSO₄ (Fisher Scientific CAS 7758-98-7), giving **1 mM final**, vortex
2. Add 63 μ L of master-mix to each sample for a total of 500 μ L reaction volume
3. Vortex to mix and allow incubation at RT for 1 hour

C. MeOH/CHCl₃ extraction to remove excess probe

1. Add 2 mL (4x volume) MeOH to 8 mL/2-dram glass vial.
2. Add 0.5 mL (1x Volume) CHCl₃
3. Add samples into the dram vial. For SILAC samples, add light samples first followed by heavy samples.

4. Add 1.0 mL* (3x volume) H₂O to sample microfuge tubes and wash out any residual proteome. Take this 1.0 mL* of H₂O and add to the dram vial already containing MeOH/CHCl₃/sample
 *for SILAC samples (1 mL total sample volume), add 1.0 mL H₂O
 *for non-SILAC samples (0.5 mL total sample volume), add 1.5 mL H₂O
5. Vortex.
6. Centrifuge 1,400 g x 3 min, 4°C; Allow dram vial to rest for 3 min or until upper layer is clear
7. Remove top and bottom layers, leaving protein interface
8. Transfer protein interface to **screw-cap microcentrifuge tube** (Sarstedt D-51588) in 600 µL MeOH
9. Rinse dram vial with 150 µL CHCl₃, and add this 150 µL CHCl₃ to the screw-cap microcentrifuge tube. Vortex.
10. Rinse dram vial with 600 µL H₂O, and add this 600 µL H₂O to the screw-cap microcentrifuge tube. Vortex
11. Centrifuge 1,400 g x 3 min
12. Remove top and bottom layers, leaving protein interface
13. Add 600 µL MeOH, sonicate with probe sonicator 3 pulses (output 3)
14. Centrifuge 14,000 rpm for 5 min to pellet protein, remove MeOH

D. Reduction, Alkylation

1. Re-suspend pellet in 500 µL of 6 M urea/25 mM ammonium bicarbonate (ambic)
*3.6g Urea (MW=60.06) in 1 mL 250 mM NH₄HCO₃ stock
 (MW NH₄HCO₃=79.06; 198 mg into 10 mL gives 250 mM stock, make fresh)
 Bring volume up to 10 mL with milliQ water*
2. Sonicate with probe sonicator 3 pulses (output 3)
3. Add 5 µL 100x (1M, 15.4 mg/100 µL in H₂O) DTT, giving 10 mM final (**make fresh**)
4. Incubate 65 °C 15 min, let cool to RT or cool briefly on ice
5. Add 40 µL 50x (0.5 M, 28 mg/300 µL in H₂O) IAA, giving 40 mM IAA (**make fresh**)
6. Incubate RT 30 min, in the dark

E. MeOH/CHCl₃ extraction to remove excess reagents from reduction/alkylation

1. Add 2 mL (4x volume) MeOH to 8 mL/2-dram glass vial.
2. Add 0.5 mL (1x Volume) CHCl₃
3. Add samples into the dram vial.
4. Add 1.5 mL (3x volume) H₂O to **screw-cap microcentrifuge tube** (Sarstedt D-51588) and wash out any residual proteome. Take this 1.5 mL of H₂O and add to the dram vial already containing MeOH/ CHCl₃/sample
5. Vortex.
 Centrifuge 1,400 x g for 3 min, 4°C; Allow dram vial to rest for 3 min or until upper layer is clear
6. Remove top and bottom layers, leaving protein interface
7. Transfer protein interface to a new **screw-cap microcentrifuge tube** (Sarstedt D-51588) in 600 µL MeOH
8. Rinse dram vial with 150 µL CHCl₃ and add this 150 µL CHCl₃ to the screw-cap microcentrifuge tube. Vortex.
9. Rinse dram vial with 600 µL H₂O and add this 600 µL H₂O to the screw-cap microcentrifuge tube. Vortex
10. Centrifuge 1,400 g x 3 min
11. Remove top and bottom layers, leaving protein interface

12. Add 600 μ L MeOH, sonicate with probe sonicator 3 pulses (output 3)
13. Centrifuge 14,000 rpm for 5 min to pellet protein, remove MeOH

F. Trypsin digestion

1. Re-suspend protein pellet in 500 μ L 25 mM ambic.
2. Sonicate with probe sonicator 3 pulses (output 3)
3. Reconstitute 20 μ g of sequence-grade trypsin (1 vial) with 40 μ L of 25 mM ambic.
4. Add 7.5 μ g (15 μ L) Trypsin/Lys-c to sample and incubate at 37°C for 3 hrs.

G. Labeled protein enrichment and elution

1. Remove 100 μ L aliquot of 50% avidin slurry per sample (in this example, 1 bed volume = 100 μ L) and place in 15 mL conical. Wash 3X by adding 10 mL of PBS, centrifuging at 1,400 x g for 1 min and decanting.
2. Resuspend washed beads in 1 bed volume of PBS. Transfer equivalent volume (approx. 130 μ L) of beads to separate 15 mL conicals. Use P1000 pipette tips to allow slurry to transfer evenly ensuring the amount of beads is equal in each tube.
3. Add the 500 μ L digested samples to the beads; washing the screw-cap tube with 1.0 mL PBS.
4. Add PBS for a total volume of 5.5 mL per sample.
5. Incubate RT, rotating for 1 hr
6. Spin 1,400 x g for 1 min to pellet; remove ~90% of supernatant by decanting **ONCE**
 ***Do not decant more than one time!!! Your beads will become re-suspended and you will dump out all of the beads on the second pour!

Wash beads/elution

1. Wash 3X with 25 mM ambic by adding 10 mL, decanting **ONCE**, and centrifuging 1,400 x g for 3 min.
2. Wash 3X with LC-MS water by adding 10 mL, decanting **ONCE**, and centrifuging 1,400 x g for 3 min.
3. Transfer beads in 200 μ L of LC-MS grade water into low-bind microcentrifuge tubes.
4. Spin 1,400 x g for 3 min to pellet, remove ~90% of supernatant (as much as you possibly can) by pipetting with gel tip.
5. Elute peptides by adding 150 μ L of elution buffer (50% ACN, 0.1% formic acid (FA)), incubating for 3 min., then spinning 1,400 x g for 3 min. Wait 1 min. to let beads settle after centrifuge has stopped. Now transfer 140 μ L (of the original 150 μ L) of eluate to a new low-bind tube using a gel loading pipette tip. Be EXTREMELY careful. DO NOT TRANSFER BEADS.
6. Repeat the elution step (Step #5) 2X, instead starting with 140 μ L of elution buffer on the 2nd and 3rd elution. Transfer all eluent to a bio-rad biospin column (to trap any beads that accidentally transferred), collect in clean microcentrifuge tube by centrifugation (1,400 x g for 3 min)
7. Poke holes from inside towards outside of microcentrifuge cap (minimize plastic debris in sample).
8. Snap freeze peptides in liquid N₂ (speeds up the drying down process) and dry down on speedvac. Do not allow liquid N₂ to enter into tube.
***AVOID performing the speedvac drying step overnight because of peptide adsorption issues described in step 10 below**
9. Re-suspend peptides in 50 μ L of mobile phase A (LC-MS water + 0.1% FA)

Stage Tip cleanup of peptides

A. Protocol to prepare the StageTip

Materials:

Empore Extraction Disks:

Octadecyl (C18)-bonded silica, part no: 2215, Vendor: 3M, St. Paul, MN


200 μ l Pipette tips

Needle: Kel-F Hub (KF), point style 3, gauge 16(1.19mm) part no: 90516

Plunger Assembly: N, RN, LT, LTN for model 1702 (25 μ l)



1. Place one Empore extraction disk into a petri dish.
2. Use the Needle to pick up the stationary phase from the disk.
3. Use the plunger to insert into the 200 μ l pipette tip.



B. C18-StageTips for clean-up and concentration

*** Use gel loading tips to be able to load the solvents and the sample right on top of the disk!**
Do not dry the disk at any time! Disk without any liquid will dry out in 1-3 min! Repeat beginning from Step 1 if you dry the disk.

1. Condition column by wetting disk with 20 μ L methanol (0.9G for 1 min)
 - a. Do not let tips dry completely
2. Add 20 μ L Buffer B (80% ACN + 0.1% FA) to stage tip. Spin at 900 x *g* for 1 min.
3. Add 20 μ L Buffer A (Water + 0.1% FA) to stage tip. Spin at 900 x *g* for 1 min.
4. Load sample directly to the disc using gel loading tip. Spin at 900 x *g* for 1 min.
5. Wash using 20 μ L Buffer A. Spin at 900 x *g* for 1 min.

- *Repeat step three times. This step might take longer than 1 minute, but check it periodically and ensure the stage tip does not dry.
6. Switch to low-bind tube.
 7. Elute using 20 μL buffer B. Spin at $900 \times g$ for 1 min.
 8. *Repeat step three times. This step might take longer than 1 minute, but check it periodically and ensure the stage tip does not dry.
 9. Spike in synthetic peptide standards for normalization and QC purposes: prepare a 100 fmol/ μL Angiotensin/Vasoactive (AV) standard solution and add 2.5 μL of peptide mix into each sample. *after re-suspending samples in 50 μL of mobile phase in step 12, the AV concentration in samples is 5 fmol/ μL and injection of 2 μL onto the LC-MS results in analysis of 10 fmol of AV peptides.
 10. Poke holes from inside towards outside of microcentrifuge cap (minimize plastic debris in sample).
 11. **Snap freeze peptides in liquid N₂** (speeds up the drying down process) and dry down on speedvac. **Do not allow liquid N₂ to enter into tube.**
***AVOID performing the speedvac drying step overnight because of peptide adsorption issues described in step 10 below**
 12. Re-suspend peptides in 50 μL of mobile phase A (LC-MS water + 0.1% FA) and analyze 3 μL of sample on LC-MS.
***note that 3 μL injection volume is generally recommended; the volume may need to be adjusted based on proteome input and protein/peptide coverage**
 13. **DO NOT store peptides in low-bind tubes because peptides will adsorb to the plastic and reduce signals on the LC-MS. After reconstituting peptides with mobile phase A, transfer 25 μL aliquots of peptides into 2X amber vials with spring-loaded glass inserts for storage. One vial is stored as backup while the other vial will allow 2X injections on the QE+ protein nanoLC-MS.**

Revised: 08/16/2022