

Preparation of HeLa peptide standards for LC-MS QC:

- 1) Take one vial of Pierce HeLa Protein digest. It should have 20 µg of tryptic peptides.
- 2) Prepare mobile phase A using 0.1% LCMS grade formic acid in LCMS grade water. Mix and test the pH on litmus paper. The color change should indicate pH ~2-3.
- 3) Make a solution of 100 ng/µL of HeLa tryptic digest. With 20 µg, one vial of lyophilized peptides should be resuspended in 200 µL of mobile phase A.
- 4) Mix everything in the vial by gently vortexing.
- 5) Aliquot 25 µL into amber vials with spring inserts. Label with the date, contents, and your initials. Store in -80 freezer.
- 6) To use HeLa peptides for QC on LCMS instruments, inject 2 µL (200 ng total) for analysis.

Revised: 07/28/2022