

Preparation of PRTC MS Standard

For 1 mL of PRTC (5 fmol/ μ L)

- 1) Add 990 μ L of water + 0.1% formic acid to a glass vial
- 2) Add 10 μ L of Pierce™ Peptide Retention Time Calibration Mixture (PRTC) (500 fmol/ μ L stock) to the vial
- 3) Store in 25 μ L aliquots in spring loaded glass inserts in amber vials @ -80 °C (please include the concentration, date prepared, and initials of lab member)

Revised: 08/11/2022